

## SPITZ AIR COMPRESSOR 230 V

END-253-750-010



### PRODUCT DESCRIPTION

Power: 0.34 kW, 21 kg max. pressure 8 bar, boiler volume 9 L

### ADDITIONAL INFORMATION

<b>Animal</b>	Companion animals
<b>SKU</b>	END-253-750-010
<b>Brand</b>	SPITZ
<b>Procedure Type</b>	Lithotripsy, Minimally invasive surgery
<b>Endoscopy Device</b>	Accessory
<b>Animal Type</b>	Companion animals

#### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

