

SPITZ STERILE FILTER WITH LL-CONNECTOR (PACKAGE OF 10 PIECES)

END-300-902-007-PK10



PRODUCT DESCRIPTION

SPITZ Sterile Filter with LL-Connector Contains 10 pieces.

ADDITIONAL INFORMATION

| | |
|-------------------------|--|
| Animal | Companion animals, Equine and Large Animal |
| SKU | END-300-902-007-PK10 |
| Brand | SPITZ |
| Procedure Type | Minimally invasive surgery |
| Endoscopy Device | Accessory |
| Animal Type | Companion animals, Equine and Large Animal |

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

