

SPITZ ADVANCED SEMI-RIGID URETERO-RENSCOPE, 425MM WORKING LENGTH

END-253-905-540



PRODUCT DESCRIPTION

SPITZ fully autoclavable Spitz Advanced Semi-Rigid Uretero-Renoscope for artificial insemination - Distal Tip: 7.0 Charr. - Distal: 8.0 Charr. - Proximal: 9.5 Charr. - Working Channel 5.8 Charr. - Field of View 65° - Direction of View 5°

ADDITIONAL INFORMATION

Animal	Mixed
SKU	END-253-905-540
Brand	SPITZ
Procedure Type	Artificial Insemination , Cystoscopy, Diagnostic Endoscopy , Uretaroscopy, Urethroscopy

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights