

**ARTHVIEW® ARTHROSCOPE 30° DIRECTION OF VIEW**

END-225-140-030

**PRODUCT DESCRIPTION****SPiTZ Fully Autoclavable Arthroscope**

- Working Length : 175mm
- Diameter: 4.0mm
- Angle of view: 30°

**ADDITIONAL INFORMATION**

<b>Animal</b>	Companion Animals,Equine, Exotics
<b>SKU</b>	END-225-140-030
<b>Brand</b>	SPiTZ

**Good to know...**

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

