

**ARTHVIEW® ARTHROSCOPE 70° DIRECTION OF VIEW**

END-225-140-070

**PRODUCT DESCRIPTION****SPiTZ fully autoclavable Arthroscope**

- Working Length : 175mm
- Diameter: 4.0mm
- Angle of view: 70°

**ADDITIONAL INFORMATION**

|               |                                   |
|---------------|-----------------------------------|
| <b>Animal</b> | Companion Animals,Equine, Exotics |
| <b>SKU</b>    | END-225-140-070                   |
| <b>Brand</b>  | SPiTZ                             |

**Good to know...**

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

