

LAPAROSCOPE 4K 30°

END-300-155-430



PRODUCT DESCRIPTION

SPiTZ fully autoclavable 4k Laparoscope - Working Length : 300mm - Diameter: 5mm - Angle of view: 30°

ADDITIONAL INFORMATION

Animal	Companion Animals,Equine, Exotics
SKU	END-300-155-430
Brand	SPiTZ
Procedure Type	Bronchoscopy, Cystoscopy, Respiratory Examination, Rhinoscopy

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.
Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!
Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.
Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights