

## LAPAROSCOPE 30°

END-300-091-130



### PRODUCT DESCRIPTION

## SPiTZ fully autoclavable Laparoscope

- Working Length : 330mm
- Diameter: 10mm
- Angle of view: 30°

### ADDITIONAL INFORMATION

|               |                                   |
|---------------|-----------------------------------|
| <b>Animal</b> | Companion Animals,Equine, Exotics |
| <b>SKU</b>    | END-300-091-130                   |
| <b>Brand</b>  | SPiTZ                             |

#### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

