

## PROVET MULTI CENTRIFUGE - CLEARANCE

146-016-CL



### PRODUCT DESCRIPTION

**Now the number 1 choice for most veterinary practices, the Pro - Vet Multi is a reliable and easy to use machine designed specifically for large animal (equine) use**

The Pro-Vet Multi has all the features of the standard Burtons Pro-Vet HE but is designed for large animal use due to its size.

Includes Swing Rotor BRK5506 & BRK5401.  
Various adaptors available but not included

#### Key Features:

- Easy to use display
- Extra thick stainless steel bowl
- Easy to clean and will stay rust free
- **Exceptional quality brushless choice of two rotors:**
- **BRK5505:** 4x100 with tube size 45x100
- 4 Programs: Sperm 1/ Sperm 2/ Faeces/ Urine
- **BRK5401:** 45x100 with tube size 2x75 & 11x40
- 5 Programs: Haematocrit/ Serum 1/ Serum 2/ Urine
- Dimensions: (HWD mm) 275x325x470
- Weight: (kg) 22 (without rotor)

Supplied With a 5401 or 5505 Combination Rotor - Holds 24 x Capillary Tube & 12x2ml Tubes

**S/N: 742428**

## ADDITIONAL INFORMATION

<b>Animal</b>	Equine
<b>SKU</b>	146-016-cl
<b>Brand</b>	Centurion

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.  
**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!  
**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.  
**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights