

WELCH ALLYN PNEUMATIC WALL DIAGNOSTIC SET

WA77764



PRODUCT DESCRIPTION

This Welch Allyn Pneumatic Wall Mounted Diagnostic Set puts everything you need in one convenient package.

Energy efficient wall organiser, with handle-based rhesostats and automatic on/off sensors, improving work flow due to contoured handles that automatically turn on and off, and eliminating cross contamination by reducing contact with no on/off switch.

The wall transformer holds 2 Welch Allyn handles and features Optisense, a reliable sensor that turns handles on when they are removed from the wall set.

Key Features:

- Help reduce cost with a more energy-efficient product
- Improve workflow with contoured handles that automatically turn on and off
- Help reduce cross-contamination by reducing contact with no switches

Kit Includes:

- Co-axial ophthalmoscope
- 3 x Specula
- Pneumatic Otoscope

ADDITIONAL INFORMATION

Animal

Small Animal

SKU	WA77764
Overall Product Dimensions (mm)	20.2cmW x 23.4cmH x 12.1cmD
Generic Size	Small
Colour/Style	White
Brand	Welch Allyn
Digital	No
Diagnostic Head Type	Pneumatic Otoscope

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights