

## MACH LED 130 PLUS EXAMINATION LIGHT

MACHLED130PLUS-BUNDLE



### PRODUCT DESCRIPTION

---

**The Mach LED 130 Plus model is a compact fixed focus examination light with powerful 100,000 lux LED technology.**

**Part of the Mach LED 130 range, the Mach LED 130 Plus is a compact fixed focus examination light which has high light intensity for critical examinations.**

Producing a cool light the LED technology is much more effective than conventional light sources, with favourable flow channels that prevent disruptive heat development. The temperature increases in the surgeon's working area is almost non-existent. Features include a faceted multi lens system, a colour rendering index of Ra=96 and superior handling design allowing the light to be positioned exactly over the wound field. The Plus also produces a brighter light for even greater illumination.

When the ceiling height is below 2.8m or above 3.0m please state the fixation height and height to underside of the false ceiling so the correct length pole can be supplied.

#### Key Features:

- Brightness (LUX): 100,000 (at a distance of 1m)
- Faceted multi lens system for homogenous and shadow-free illumination
- Colour Rendering Index (CRI): Ra = 96 for accurate tissue colour identification
- Authentic colour rendering to reduce eye strain
- Superior handling design - allowing light to be positioned where needed
- 0.5°C temperature increase
- Colour temperature (Kelvin): 4500
- Dimming range (%): 50-100
- Number of LED: 19
- Size of the light field (cm): 12 (fix) - Focusable version see (130F model)
- Total power consumption (W): 16
- Lifespan of the LED's (h): 60.000
- Working distance (in cm): 70-140

- Chassis of the light (cm): 33

### Hygienic

High quality materials provide security and avoid infiltrations of dust, dirt, and liquids inside the lamp head. The disk sealings of the light outlets and the circumferential sealing cord guarantee a perfect product down to the smallest detail.

### Effective

Powerful LED-technology helps keeping a cool head through an almost non-existent heat development.

## ADDITIONAL INFORMATION

|                               |  |
|-------------------------------|--|
| <b>Animal</b>                 | N/A  |
| <b>SKU</b>                    | MACHLED130PLUS-bundle                        |
| <b>Colour/Style</b>           | White  |
| <b>Brand</b>                  | Dr Mach                                      |
| <b>Autoclavable Handle</b>    | Yes  |
| <b>Colour Temperature (k)</b> | 4500   |
| <b>Focusable</b>              | No   |
| <b>Head Type</b>              | Single                                       |
| <b>LED Life (hours)</b>       | 60,000                                       |
| <b>Light Source</b>           | LED  |
| <b>Lux @1 meter</b>           | 100000                                       |
| <b>Produceure</b>             | Dental, Exam, Minor-Operation, Task, Theatre |

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.  
**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!  
**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.  
**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights