

## HUGEMED BR-M40H BRONCHOSCOPE - FOR FELINE/CANINE/ZOO - 600MMWL X 4MMOD WITH 1.2MM CHANNEL

END-BR-M40H



### PRODUCT DESCRIPTION

**The BR-M40H Semi-Disposable Flexible Bronchoscope is a state-of-the-art tool engineered for the diagnosis and treatment of respiratory conditions in veterinary patients.**

**Patient Type:** Feline / Canine / Zoo

**Procedure Type:** Upper respiratory - Pathological samples / biopsies / removal of foreign bodies / BAL / Cytology

Semi-disposable scopes are a cost-effective solution to limiting contamination risks, incurring less processing costs, and reducing time spent on maintenance whilst remaining the ideal solution for procedures.

Combining advanced design with high-quality imaging, it is perfectly suited to veterinary use, delivering exceptional precision in airway examinations, biopsies, foreign body retrieval treatments, and surgeries. This versatile bronchoscope is specifically designed for procedures involving the nasal cavity, upper respiratory tract, trachea, and bronchus, ensuring reliable performance across a range of diagnostic and therapeutic applications.

It has a working length of 600mm with 2-way rotation making it the ideal choice for veterinarians. The insertion tube outer diameter for the endoscope is 4mm, with a minimum working channel inner diameter of 1.2mm. The imaging quality is at Full HD resolution (1280 x 720 pixels) for crystal-clear visuals to help you navigate challenging procedures. The field of view for the scope is  $120^{\circ} \pm 10\%$ , with a  $0^{\circ} \pm 10^{\circ}$  forward-facing direction of view and a depth of field ranging from 3mm to 50mm.

The BR-M40H provides high-resolution imaging via the MS-8 Video Processor, ensuring clear, detailed visuals for examinations of the nasal cavity, upper respiratory tract, trachea, and bronchus. The unit can be used on small animal, wildlife, zoo, avian and exotics. It can be utilised in a mixed practice or even in a remote setting.

The lightweight handle of the unit includes a suction port for seamless integration with irrigation systems and a biopsy port for

channel instruments, with 2 customisable buttons that can be programmed from the all-in-one MS-8 processor to control photographing, videoing, white balance adjustments and graphics freezing with the selected buttons displayed as hotkeys.

There is a manual lever at the base of the handle which controls the bending section and the distal tip of the scope for precise navigation. The high-quality medical grade patient tube is equipped with an integrated camera and LED light source. These components deliver high-definition visuals, transmitting image data to the MS-8 Image Processor via a 2500mm connecting wire.

Every bronchoscope is packed in sterilised, single-use packaging with the date of manufacture and expiration date clearly indicated on the label, ensuring an aseptic shelf life of 3 years. This minimises cross-contamination risks, enhancing safety for both patient and practitioner. Professional disinfectant wipes that are effective against bacteria, yeasts, viruses, and fungi which are tested to veterinary standards (European PT3 test standards for microbiological testing on key veterinary equipment) are suitable for cleaning the outside surfaces of these scopes.

With seven models within the BR series, there is ample choice to meet the unique and diverse needs of veterinary practitioners who are making diagnosis and undertaking treatments on patients. The various sizes satisfy the requirements of different operational settings. This range is best suited to alveolar lavage, airway examinations, biopsies, foreign body removal, intubation, medication delivery and secretion suction procedures. The units are fast to set up, with the ability to run on battery for over 4 hours, making them the ideal choice for out of hours and remote work.

The BR-M40H bronchoscope offers reliability, flexibility, and high performance, making it a trusted choice for veterinary professionals seeking excellence in endoscopic care. It is compatible with all the conventional accessories and 1.0mm instruments used in endoscopic care.

### Key Features

- Full HD Resolution for clear, high-quality imaging.
- Built-in LED light source and high-definition camera.
- Variable field of view and sharp focus in narrow spaces and at great depths.
- High-quality patient tube with flexible distal end for smooth navigation.
- Supports a wide range of procedures versatile applications and advanced techniques.
- Easy and wide aspiration and sampling through the working channel.
- User-friendly multi-functional handle with programmable hotkeys for ease of use.
- Manual lever for precise control and adjustment when using the bending section.
- Safe and hygienic sterilised packaging with a 3-year aseptic shelf life, reducing cross-contamination risks.

### Specification

<b>Endoscope Model</b>	BR-M40H - Semi-Disposable Flexible Video Bronchoscope
<b>Field of View (°)</b>	120° ±10%
<b>Direction of View (°)</b>	0°±10°
<b>Depth of Field (mm)</b>	3mm - 50mm
<b>Insertion tube diameter(mm)</b>	4.0mm
<b>Working Length (mm)</b>	600mm±10%
<b>Working Channel (mm)</b>	1.2mm±10%
<b>Bendable Portion (°)</b>	210° (upwards) and 210° (downwards), with a deviation of -5% (excluding upper limit)
<b>Functions</b>	Photographing / Videoing / White Balance / Graphics Freezing
<b>Connections</b>	USB / HDMI
<b>U/D, L/R rotation(ways)</b>	2 Way Rotation
<b>CMOS Pixel (px)</b>	1280px*720px
<b>Optical</b>	Full HD, Wide Angle

**Overall Length (mm)** 860mm±10%

**Operational Environment Temperature (°C)** 5-40°C(41-104°F)

**Supporting Products:**

- END-310-901-120 - Biopsy forceps Long Oval Cup, fenestrated D: 1.0 mm, 120 cm length
- END-310-960-120 - Alligator Forceps flex. foreign body, Forceps D: 1 mm, WL 120 cm
- END-NB3FR120 - 1.0mm x 1200mm Basket
- END-SG3FR120 - 1.0mm x 1200mm 3 Prong Grasper
- END-RC3FR120-S - 1.0mm x 1200mm pack of 10 sterile BAL catheters
- END-253-962-003 - Brush for Cytology, 3 Charr., For single use
- END-310-901-120 - Biopsy forceps Long Oval Cup, fenestrated D: 1.0 mm, 120 cm length
- END-17021182 - Mucus Collection Trap 30ml - Pk of 10
- END-207230420 - Cleaning Brush Ø2mm, length 2300mm, Disposable Double-Sided, Ball Head, Tube Ø1.0mm

**List of Other Semi-Disposable BR Scopes**

Semi-disposable BR scopes							
Model	BR-M58H	BR-M58H-1500	BR-M52H	BR-M52H-1000	<b>BR-M40H</b>	BR-M32	BR-M25
<b>Insertion tube diameter(mm)</b>	6.0	6.0	5.0	5.0	<b>4.0</b>	3.1	2.5
<b>Working channel(mm)</b>	2.8	2.8	2.2	2.2	<b>1.2</b>	1.2	1.2
<b>Working length(mm)</b>	600	1500	600	1000	<b>600</b>	600	600
<b>U/D, L/R rotation(ways)</b>	2	4	2	4	<b>2</b>	2	2
<b>CMOS Pixel</b>	1280*720 HD	1280*720 HD	1280*720 HD	1280*720 HD	<b>1280*720 HD</b>	400*400	400*400

**ADDITIONAL INFORMATION**

<b>Animal</b>	Companion animals
<b>SKU</b>	END-BR-M40H
<b>Overall Product Dimensions (mm)</b>	860
<b>Colour/Style</b>	White
<b>Brand</b>	HugeMed
<b>Procedure Type</b>	Bronchoscopy, Diagnostic Endoscopy , Respiratory Examination, Rhinoscopy, Endoscopy, Laparoscopy, Arthroscopy, Cystoscopy, Tracheoscopy, Bronchoscopy, Cystoscopy, Otoscopy
<b>Biopsy Channel</b>	1.2mm
<b>Working Length (mm)</b>	600
<b>Fitting</b>	Universal

## Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights