

**ACEM ST1 EXAM LIGHT - 80CM POLE - BLUE - CLEARANCE**

ACEMST1.3/RC-80-BLUE-CL

**PRODUCT DESCRIPTION**

ACEMST1 Ceiling Light - 80cm - Blue

ACEMST1.3 Arm Ceiling Examination Light  
For Veterinary Use-Standard Ceiling Pole 80cm  
For floor to ceiling height: 270cm-280cm  
Colour - Blue  
Articulated Arm  
100,000 Lux @ 50cm  
Colour Temp: 3,000K to 5,500K  
Light Field Diameter: 26cm  
LED Average Life: 60,000h (About)

**ADDITIONAL INFORMATION**

<b>Animal</b>	N/A
<b>SKU</b>	ACEMST1.3/RC-80-BLUE-cl
<b>Colour/Style</b>	Blue
<b>Brand</b>	Acem
<b>Colour Temperature (k)</b>	3,000K to 5,500K
<b>LED Life (hours)</b>	60,000h
<b>Lux @1 meter</b>	100,000 Lux @ 50cm

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights