

MACH LED CEILING LIGHT - FOCUSABLE - STERILISABLE HANDLE - 25CM - CLEARANCE

MACHLED130F-C-SH-25-CL

**PRODUCT DESCRIPTION**

MACH LED Ceiling Light - Focusable - Sterilisable Handle - 25cm

MACH LED 130F Ceiling Mounted Light
with 25cm Ceiling Pole Kit
Includes Removeable Sterilizable Handle

Please note: Anchorage Ring is not included

ADDITIONAL INFORMATION

Animal	N/A
SKU	MACHLED130F-C-SH-25-cl
Brand	Dr Mach
Autoclavable Handle	Yes
Focusable	Yes

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights