

BESPOKE LO-LIFT TABLE TOP (170CM L X 80CM W) - CLEARANCE

XX-2057



PRODUCT DESCRIPTION

Bespoke Lo-Lift Table Top
Bespoke Length (170L x 80 W)

The lo-lift table series built by Burtons has been designed around the Vet, with ultra low minimum height, 3 pedal braking , electric battery option and high quality frame. This table is a staple of many veterinary clinics. With a bespoke table top fitting this table can be configured for your needs allowing for most of the Burtons table top range to be fitted with additional heated and lead lined options available . Simply select the top you need to fit you Lo-Lift Base.

ADDITIONAL INFORMATION

Animal	Mixed
SKU	XX-2057
Brand	Burtons

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights