

22MM INTERSURGICAL - FLEXTUBE™ CIRCLE BREATHING SYSTEM

103-064



PRODUCT DESCRIPTION

This circle breathing system consists of 2 limbs, 22mm in diameter with swivel Y-piece connectors, a luer lock port for the monitoring line, and a safety cap. This ensures optimal positioning for accurate and reliable gas monitoring, minimising the possibility of the monitoring line becoming caught or entangled.

Along with this, the limbs corrugated tubing also contains a 22M straight connector to convert one of the limbs from a 22F to a 22M connection for an airtight fit, minimising the risk of disconnection or occlusion.

This product will fit the Cyclo-Flo and the Semi Disposable Circle Circuit.

- 22mm bore with 22mm cuffs on each tube
- 1.6 metre length
- Swivel Y-piece patient connector (22M/15F)
- Made in the UK to UKAS and ISO standards

ADDITIONAL INFORMATION

Animal	Small Animal
SKU	103-064
Tubing Length	1.6 Meter
MRI Compatible	No
Circuit Type	Parallel

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights