

## WELCH ALLYN HANDLES

WA71020-B-1



## PRODUCT DESCRIPTION

### High quality, long lasting Diagnostic handles

#### Choice of:

- Standard 3.5v rechargeable handle
- Non Rechargeable handpiece for C cells
- 3.5v Rechargeable convertible handle (normal batteries)
- Lithium Ion Smart Handle - Powers any instrument head for 110 minutes - thats twice as long as a conventional handle

#### Key Features:

- Sure grip smooth, non-abrasive knurled finish ensures a no-slip grip
- Chrome-plated brass handles with durable, long-life construction
- Safety switch to prevent accidental draining of battery
- Metal head connectors for durability

## ADDITIONAL INFORMATION

Animal

N/A

<b>SKU</b>	WA71020-B-1
<b>Generic Size</b>	X-Small
<b>Digital</b>	No
<b>Power Type</b>	3.5v Rechargeable, Lithium Rechargeable, Non-Rechargeable

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.  
**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!  
**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.  
**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights