

WELCH ALLYN HANDLES

WA71020-B-1



Price: £104.96 - £185.00
Date Printed: 28/03/24

Category: Diagnostic Imaging > Diagnostic Essentials > Diagnostics

PRODUCT DESCRIPTION

High quality, long lasting Diagnostic handles

Choice of:

- Standard 3.5v rechargeable handle
- Non Rechargeable handpiece for C cells
- 3.5v Rechargeable convertible handle (normal batteries)
- Lithium Ion Smart Handle - Powers any instrument head for 110 minutes - thats twice as long as a conventional handle

Key Features:

- Sure grip smooth, non-abrasive knurled finish ensures a no-slip grip
- Chrome-plated brass handles with durable, long-life construction
- Safety switch to prevent accidental draining of battery
- Metal head connectors for durability

ADDITIONAL INFORMATION

| | |
|-------------------|---|
| Animal | N/A |
| SKU | WA71020-B-1 |
| Size | X-Small |
| Digital | No |
| Power Type | 3.5v Rechargeable, Lithium Rechargeable, Non-Rechargeable |

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.
Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!
Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.
Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights