

## OPTICLAR LED DIAGNOSTIC WALL STATION - WITH SLIT OTOSCOPE & OPHTHALMOSCOPE

OPT700.041.010



### PRODUCT DESCRIPTION

---

## Veterinary Diagnostic Set - Wall Mounted - Mains & C-Cell

Opticlar's VisionVet range is designed to simplify and enhance diagnostic examinations across species, from exotics to companion animals and larger breeds.

Engineered for durability and performance, each unit features TrueTone LED modules, delivering exceptional illumination with a long lifespan. The 10-year LED module guarantee means maintenance costs are significantly reduced over the lifespan of the equipment. The high-quality optics ensure excellent visualisation, while the robust construction offers lasting reliability.

An automatic on/off function and compatibility with traditional manufacturer connections help protect your existing equipment investments. Additionally, the multi-voltage design allows seamless switching between cell power and rechargeable options without the need to replace bulbs.

### Set includes:

- S1 Veterinary LED slit otoscope
- L28 Practice ophthalmoscope with cobalt blue filter, light sensitive filter and 5 apertures
- Twin handle wall unit with auto on/off and 12ft approx cord handles
- Twin port charger
- 3 metal tips - 57.0mm x 4.0mm diameter, 65.0mm x 6.0mm diameter, 90.0mm x 7.0mm diameter with German brand compatible fitting

### Key Features

- TrueTone multi-voltage LED slit head
- Can be 'locked' on to handle
- Excellent LED illumination
- Bayonet tip connection
- Swivelling 3x magnification lens
- Durable, solid metal construction
- Opticlar and Welch Allyn connection
- 10 year LED warranty
- 5 year product warranty
- Supplied with 3 metal tips

### Veterinary LED slit Otoscope

- TrueTone multi-voltage LED illumination
- Distal fibre optics gives cool light to examination site with no reflections

### L28 Practice LED Ophthalmoscope

- Polarising filter practically eliminates corneal reflex and can be used on all apertures
- Sealed optics for long life and reduced running cost
- 28 lenses from -25 to +40d
- Independently selectable red free filter for use on all apertures provides increased contrast
- 18 filter/aperture combinations for increased diagnostic capabilities
- Cobalt blue filter to help detect corneal abrasions
- Easy to operate lens wheel for accurate focusing
- Smooth action thumb wheel for aperture selection

### Wall Station

- Designed without flat surfaces that can attract unwanted items such as used tips etc.
- Rheostat brightness control on handle
- Handles have Auto on/off facility
- Spiral cords have a reach of approx. 3mtr/ 10ft
- Locking collar secures instrument head in to cord handles guarding against loss and theft
- Wall units can be used as a direct replacement for other leading makes either directly or by use of simple adaptors
- Can be supplied with mobile trolley unit for maximum portability

### Specifications:

- LED illumination life: 10,000 hours
- LED guarantee: 10 years

### ADDITIONAL INFORMATION

<b>Animal</b>	Mixed
<b>SKU</b>	OPT700.041.010
<b>Brand</b>	Opticlar
<b>Diagnostic Head Type</b>	Co-axial Ophthalmoscope, Slit Otoscope

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights