

OPTICLAR LARYNGOSCOPE MILLER BLADES - WITH MAGNIFICATION - SIZE 0 - 55 X 10MM

OPT600.020.000VMAG



PRODUCT DESCRIPTION

Veterinary Laryngoscope Miller Blades - With Magnification - Size 0 - 55 x 10mm

Opticlar have developed a unique range of laryngoscopes and other ancillary intubation accessories designed to aid and simplify airway management in exotics, companion animals and larger species. Blades only

Key Features

- Miller Vet Blades with 3x Magnification
- Stainless construction, easy to clean and fully autoclavable
- Size 0
- 55 x 10mm
- Fibre optic light transmission prevents heat build up at blade end
- Reliable and long lived with 10 year product warranty
- Compatible with ISO 7376 Green laryngoscope specification

ADDITIONAL INFORMATION

Animal	Mixed
SKU	OPT600.020.000VMag
Brand	Opticlar

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights