

**OPTICLAR METAL SPECULA - FELINE - 32MM X 3.0MM**

OPT700.000.000-3S

**PRODUCT DESCRIPTION****Veterinary Metal Specula - Feline - 32mm x 3.0mm**

Metal specula with secure double bayonet fitting for OPTICLAR VET instruments and similarly-equipped VET otoscopes from some other manufacturers.

**Key Features**

- Compatible with Pro-VET LED Otoscopes and LED Pocket Otoscopes
- Autoclavable
- Metal construction, slit or closed designs
- Bite-resistant

**ADDITIONAL INFORMATION**

<b>Animal</b>	Mixed
<b>SKU</b>	OPT700.000.000-3S
<b>Brand</b>	Opticlar

## Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights