

OPTICLAR METAL SPECULA - STENOTIC EAR - 70MM X 5.0MM

OPT700.000.000-5L

**PRODUCT DESCRIPTION****Veterinary Metal Specula - Stenotic Ear - 70mm x 5.0mm****For use with Veterinary Slit & Pocket Otoscopes**

Metal specula with secure double bayonet fitting for OPTICLAR VET instruments and similarly-equipped VET otoscopes from some other manufacturers. This ear tip has been designed to use when examining stenotic ear canals. Its length and smaller diameter make it significantly easier when introducing into a stenotic ear canal. Despite its small diameter this tip provides very good visualisation.

Key Features

- Compatible with Pro-VET LED Otoscopes and LED Pocket Otoscopes
- Autoclavable
- Metal construction, slit or closed designs
- Bite-resistant

ADDITIONAL INFORMATION

Animal	Mixed
SKU	OPT700.000.000-5L
Brand	Opticlar

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights