

## OPTICLAR FELINE & SMALL ANIMAL METAL SPECULA

700.000.000-3S-GROUPED



### PRODUCT DESCRIPTION

---

## Feline and Small Animal Specula

#### For use with Veterinary Slit & Pocket Oscopes

- Opticlar 32.0mm x 3.0mm Feline Tip (700.000.000-3S)
- Opticlar 32.0mm x 4.0mm Feline Tip (700.000.000-4S)
- Opticlar 32.0mm x 5.0mm Feline Tip (700.000.000-5S)

Feline and Small Animal Specula - These reduced length and diameter tips are ideal for use with feline and other small animals. These have a smooth exterior finish and a matt internal finish that reduces glare and maximises the amount of light thrown from the tip end. Fully autoclavable.

#### Specifications:

- Metal construction
- Autoclaveable
- Bite-resistant

### ADDITIONAL INFORMATION

SKU

700.000.000-3S-grouped

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights