

## MOBILE INFUSION STAND

I4110017



## PRODUCT DESCRIPTION

### Easy to assemble mobile infusion (drip) stand\*

Part of the Burtons value range, this practice essential with its chrome finish looks great and is easy moved around the practice thanks to its lightweight construction and 5 star castor base.

- Height adjustable (135-215 cm)
- Manoeuvrable and lightweight
- Chrome finish
- 4 Hook design
- 5 star castor base

**\*Infusion Pump Not Included**

## ADDITIONAL INFORMATION

<b>Animal</b>	Mixed
<b>SKU</b>	I4110017
<b>Colour/Style</b>	Silver
<b>Brand</b>	Burtons

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights