

## MINDRAY MR401B RE-USABLE TEMP PROBE, ADULT,

UMEC-125-000143-00



### PRODUCT DESCRIPTION

Mindray MR401B Re-usable Temp Probe, Adult,. Change when broken. Compatible with Mindray uMEC & ePM12 monitors.

### ADDITIONAL INFORMATION

<b>Animal</b>	N/A
<b>SKU</b>	UMEC-125-000143-00
<b>Brand</b>	Mindray
<b>Accessory Parameter</b>	Temp
<b>Monitor Compatibility</b>	Mindray ePM12, Mindray uMEC
<b>Change Frequency</b>	When Broken

#### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.  
**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!  
**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.  
**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

