

## SUPER STRENGTH GIVING SET

129-015\*BUNDLE



### PRODUCT DESCRIPTION

## Burtons super strength giving set with anti-kink design

Our improved giving set is now even more versatile and stronger! We've added an anti-kink section nearest the patient to further prevent occlusions and a new 'C' Clip flow shut off clamp to temporarily shut off the flow. The set is ideally suited for use with the IVY, Niki IV or VetPro and other infusion pumps or as a gravity fed line.

**Why not get your infusion pump calibrated by our workshop and ensure accurate readings, call us to discuss!**

- Overall length: 2.4 meters .O.D. 4.1mm
- 2 x 'Y' style injection ports
- 1 x Precision V Roller
- Rotating luer lock and slip luer patient connector
- 20 drop/ml dropper
- Super strong tube

### ADDITIONAL INFORMATION

<b>Animal</b>	Mixed
<b>SKU</b>	129-015*bundle
<b>Brand</b>	Burtons
<b>Catagory Sub Heading</b>	Bulk Buy Packs Available

## Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights