

DRAMINSKI LOW FREQUENCY CONVEX PROBE

DRAM-CON-LOW



PRODUCT DESCRIPTION

For pregnancy diagnostics in sheep and goats. Abdomen and thorax of cattle and horses. Horse musculoskeletal system. Umbilical region in calves.

Specifications:

- Frequency range: **2.0 - 5.0 MHz**
- Central frequency: **3.5 MHz**
- Scanning depth: **60 - 350 mm**
- Lens radius: **60 mm**
- Number of crystals: **128**

ADDITIONAL INFORMATION

Animal	Equine and large animal
SKU	DRAM-CON-LOW
ultrasound probes compatibility	Draminski BLUE
Ultrasound Probe Uses	Pregnancy diagnostics in sheep and goats. Abdomen and Thorax of cattle and horses. Horse musculoskeletal system. Umbilical region in calves

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights