

DRAMINSKI ABDOMINAL HIGH FREQUENCY LINEAR PROBE

DRAM-LINEAR-HIGH



PRODUCT DESCRIPTION

For tendon examination in a horse or lung and abdomen examination of small animals

Specifications

- Frequency range: 6.0 – 14.0 MHz
- Central frequency: 10.0 MHz
- Scanning depth: 40 – 100 mm
- Probe face length: 40 mm
- Number of crystals: 128

ADDITIONAL INFORMATION

Animal	Companion animals, Equine and Large Animal
SKU	DRAM-LINEAR-HIGH
ultrasound probes compatibility	Draminski BLUE
Ultrasound Probe Uses	Equine tendon examinations, Small animal lung and abdomen examinations

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights