

BLT E30 VET ECG MACHINE

E30VETECG



PRODUCT DESCRIPTION

The E30 Vet ECG Machine is designed to make cardiac monitoring simple, accurate, and reliable in any veterinary practice.

With a lithium-ion battery and intuitive touchscreen interface, it's easy to use whether you're in the clinic or on the move.

Key Features

- Touchscreen display for quick and intuitive operation
- 7-lead ECG cable included
- Weighs only 1.6kg, so it can be easily moved around your practice
- Built-in high-resolution thermal printer for immediate results
- Rechargeable lithium-ion battery with up to 4 hours of continuous use
- Advanced ECG algorithm with automatic measurement and baseline adjustment for fast, accurate readings
- Powerful digital filter minimises interference to produce clean, precise ECG waveforms
- USB port for efficient data export and sharing
- Cat, dog and equine pre-sets already set up
- Includes 1 x roll of printer paper, probes, clips and 7-lead ECG cable

Specifications

Size: 315mm x 215mm x 75mm
Net weight: 1.6kg
Display: 4:3" TFT LCD, resolution 480x272
Power: AC 100-240V, 50/60Hz or DC rechargeable lithium-ion battery (11.1V, 2Ah)

Safety

Class I, Type CF product with defibrillation protection
Isolation voltage: 4000V
Polarisation resistance: $\pm 500\text{mV}$
Recording
Thermal printing system with 80mm x 20m reel paper
Resolution: Vertical 8 dots/mm, Horizontal 40 dots/mm
Recording speeds: 5-50mm/s

Heart Rate

Measurement range: 30-300 bpm
Accuracy: ± 1 bpm (average 10s)

Storage

Holds up to 200 ECG data groups in internal memory
External storage supported via USB

This reliable ECG unit ensures quick, accurate cardiac diagnostics with user-friendly operation—an essential tool for busy veterinary teams looking for dependable results every time.

ADDITIONAL INFORMATION

Animal	Mixed
SKU	E30VETECG
Screen Type	Touchscreen
Power Source	Rechargeable

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.
Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!
Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.
Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights