

STAINLESS STEEL WATER DISTILLER

139-100



PRODUCT DESCRIPTION

Produce distilled water in practice with this compact stainless steel water distiller

Distilled water is free from the minerals and impurities that cause scale, blockages, and corrosion, making it ideal for autoclaves, dental units, and other equipment that rely on clean steam and water. With a one-gallon output in approximately 5.5 hours, it helps maintain consistent sterilisation, reduces avoidable downtime, and supports the long-term protection of your equipment.

Key features:

- Stainless steel boil chamber, dome, and condensing coils for robust, long-term use in clinical environments
- Produces distilled water suitable for a wide range of practice applications, from steriliser reservoirs to instrument rinse water
- Integrated safety thermostat automatically switches the unit off at the end of the cycle
- 1 gallon (approx. 3.8 litres) output every 5.5 hours, supporting routine daily demand in smaller practices or satellite sites.
- Dimensions: 20 cm wide x 36 cm tall (excluding collection bottle), weighing just 3.5 kg

What's included:

- Stainless steel water distiller unit
- BPA-free polypropylene collection bottle.
- Detachable power cord and filter housing.
- 250g distiller residue cleaner to support straightforward ongoing maintenance
- 1-year warranty

ADDITIONAL INFORMATION

SKU

139-100

Animal

N/A

Overall Product Dimensions (mm)

20cmW x 36cmH

Generic Size

Small

Colour/Style

Stainless Steel

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights