

MACH LED 115 MOBILE EXAMINATION SPOT LIGHT

MACHLED115-M



PRODUCT DESCRIPTION

The Mach Mobile examination spot light is sturdy, flexible and easy to move. Ideal for minor procedures and small rooms.

- Ultra compact design, ideal for small rooms.
- LED technology has a minimum life-span of 40,000 hours
- Easy to manoeuvre flexible gooseneck arm
- Authentic colour rendering to reduce eye strain (rendering index R(a) 95)
- Desk clamp option
- 60,000 lux at 0.5m distance
- Light head: 22cm
- Number of LEDs: 7
- Total power consumption: 10 watts

ADDITIONAL INFORMATION

Animal	N/A
SKU	MACHLED115-M
Colour/Style	White
Brand	Dr Mach
Autoclaveable Handle	No
Battery Backup	No
Colour Control	No
Focusable	No
Head Type	Single
LED Life (hours)	40,000

Light Source

LED

Lux @1 meter

30,000

Producedure

Exam, Task

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.
Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!
Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.
Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights