

MINDRAY 3-LEAD ECG WIRES, ADULT / PAEDIATRIC, SNAP, TPU, IEC (1M)

IMEC-0010-30-42733



PRODUCT DESCRIPTION

Mindray 3-Lead ECG Wires, Adu/Ped, Snap, TPU, IEC (1m). Change when broken. Compatible with Mindray uMEC & Mindray ePM12 monitors.

ADDITIONAL INFORMATION

Animal	N/A
SKU	IMEC-0010-30-42733
Brand	Mindray
Accessory Parameter	ECG
Monitor Compatibility	Mindray ePM12, Mindray uMEC
Change Frequency	When Broken

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.
Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!
Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.
Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

