

PETMAP BP EXTENSION HOSE FOR 6.5 - 13.0CM CUFF

129-578



PRODUCT DESCRIPTION

PetMap BP Extension Hose for 6.5 - 13.0cm Cuff. Change when broken. Compatible with PetMap G2 & PetMap G3 monitors.

ADDITIONAL INFORMATION

| | |
|------------------------------|--------------------|
| Animal | N/A |
| SKU | 129-578 |
| Accessory Parameter | NIBP |
| Monitor Compatibility | Vetronic Lightning |
| Change Frequency | When Broken |

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights