

## MACH LED 130 MOBILE EXAMINATION LIGHT

MACHLED130-M-GROUP\*



### PRODUCT DESCRIPTION

#### Mach LED 130 Mobile Examination Light with UK Plug

The Mach LED 130 is a compact veterinary examination light with LED technology.

Producing a cool light the LED technology is much more effective than conventional light sources. The temperature increase in the surgeons working area is almost non-existent. Features include a faceted multi-lens system, a colour rendering index of Ra=95 and superior handling design allowing the light to be positioned exactly to the wound field.

Data Common To All Models		
Technical Specifications	Mach LED 130 F	Mach LED 130
Light Intensity at 1 Meter Distance (lux)	60,000	50,000
Colour Rendering Index R(a) General	95	95
Light-Field Size (cm)	14 - 25 Focusable	Fixed Focus 17
Colour Temperature	4,500	4,500
Electronic Light Intensity Control On Lamp Head	50 - 100%	50 - 100%
Power Consumption	28 W	28 W
Number of LED's	19	19
Mains	230 V 50 Hz	230 V 50 Hz
Working Distance (cm)	70 - 140	70 - 140
Diameter of Light Head	33	33
Height Adjustment (cm)	123	123

Weight (Ceiling, Wall, Mobile)

15 Kg, 8 Kg, 20Kg

## ADDITIONAL INFORMATION

<b>Animal</b>	Small Animal
<b>SKU</b>	MACHLED130-M-group*
<b>Focusable</b>	No
<b>Head Type</b>	Single
<b>Light Source</b>	LED
<b>Lux @1 meter</b>	60,000
<b>Produceure</b>	Dental, Exam, Minor-Operation, Task

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.  
**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!  
**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.  
**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights