

BURTONS VERTICAL INFUSION POLE FOR ACRYLIC CAGE DOORS

165-490-GROUPED*



PRODUCT DESCRIPTION

A practical and efficient upgrade to your cage banks

The Burtons Vertical Infusion Rail is designed specifically for cages with acrylic doors. By positioning infusion pumps at a more accessible working height, rather than above the cages on traditional socket rails, this system enhances usability and day-to-day workflow.

The Vertical Infusion Rail provides a smart, ergonomic solution for both new acrylic door installations and can be added to your existing cages. Only compatible with Burtons Cages.

ADDITIONAL INFORMATION

SKU

165-490-grouped*

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. **Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! **Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. **Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

