

MACH LED 130 + WALL MOUNTED EXAMINATION LIGHT

MACHLED130PLUS-W



PRODUCT DESCRIPTION

Mach LED 130 + Wall Mounted Examination Light with UK Plug

The Mach LED 130 + W is a compact, wall mounted examination light with LED technology and has an extra 40,000 lux than the standard Mach 130. It produces a cool light and the temperature increase in the surgeons working area is almost non-existent.

Key Features:

- 100,000 Lux (at a distance of 1m)
- Facetted multi lens system
- Colour rendering index of of Ra=95
- Authentic colour rendering to reduce eye strain
- Superior handling design - allowing light to be positioned where needed
- 30,000 lux at 0.5m distance
- 0.5 temperature increase
- 40,000 hours bulb life
- Fixed focus

Facetted Multi-lens System

A multitude of computer-calculated facetted lenses guarantees homogeneity and lowest shadiness in the light field. Separately arranged optical systems, each with one LED module, generate their own light field, which increases the contrast effect of the light. Light intensities of 50.000 Lux can be attained without difficulty.

Superior Colour Rendition

With outstanding colour rendering indexes Ra = 95 the surgeon recognizes clearly the tiniest nuances of colour in tissue. The colour spectrum of the wound is rendered naturally with richcontrast. The OT-light clearly provides welcome relief for your eyes.

✖ Keypad on the Lamp Housing

The following light functions can be adjusted electronically, such as:

- Switching ON and OFF
- Electronic light intensity control

ADDITIONAL INFORMATION

Animal	N/A
SKU	MACHLED130PLUS-W
Colour/Style	White
Brand	Dr Mach
Autoclaveable Handle	Yes
Battery Backup	No
Colour Control	No
Colour Temperature (k)	4300
Focusable	No
Head Type	Single
LED Life (hours)	50,000
Light Source	LED
Lux @1 meter	100,000
Procedure	Dental, Exam, Minor-Operation

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights