

RESPIVET RESPIRATORY MONITOR

117-950



PRODUCT DESCRIPTION

An easy to use and low-cost monitoring solution

The RespiVet respiratory monitor unit has variable sensitivity settings and adjustable audible and visual apnoea alarms.

Features

- Apnoea Alarm settings: 10-60 seconds
- Adjustable sensitivity for patients: 0.5 to 1,000 kg
- Photo eye for handsfree alarm reset
- LED readout of respiration rate
- Acoustic signal for each breath
- Connects endotracheal tube to the anaesthetic circuit
- Patient Hose Connector: Outer : 22 mm, Inner: 15 mm
- Tracheal Tube Connector: Outer : 15 mm
- Rechargeable battery: Up to 70 hours of use
- Dimensions (cm): 13 x 13 x 6 (L x W x H)

ADDITIONAL INFORMATION

| | |
|---------------------------|---------------------|
| Animal | Mixed |
| SKU | 117-950 |
| Colour/Style | White |
| Monitor Parameters | Respiration |
| Screen Type | Digital Display |
| Power Source | Mains, Rechargeable |
| Optional Modules | No |

SD Card

No

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights