

## RESPIVET RESPIRATORY MONITOR

117-950



### PRODUCT DESCRIPTION

## An easy to use and low-cost monitoring solution

The RespiVet respiratory monitor unit has variable sensitivity settings and adjustable audible and visual apnoea alarms.

#### Features

- Apnoea Alarm settings: 10-60 seconds
- Adjustable sensitivity for patients: 0.5 to 1,000 kg
- Photo eye for handsfree alarm reset
- LED readout of respiration rate
- Acoustic signal for each breath
- Connects endotracheal tube to the anaesthetic circuit
- Patient Hose Connector: Outer : 22 mm, Inner: 15 mm
- Tracheal Tube Connector: Outer : 15 mm
- Rechargeable battery: Up to 70 hours of use
- Dimensions (cm): 13 x 13 x 6 (L x W x H)

### ADDITIONAL INFORMATION

<b>Animal</b>	Mixed
<b>SKU</b>	117-950
<b>Colour/Style</b>	White
<b>Monitor Parameters</b>	Respiration
<b>Screen Type</b>	Digital Display
<b>Power Source</b>	Mains, Rechargeable
<b>Optional Modules</b>	No

SD Card

No

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.  
**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!  
**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.  
**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights