

## HANDI+ OXYGEN ANALYSER

129-555



### PRODUCT DESCRIPTION

## So you have an oxygen therapy cage, and the oxygen concentration you provide is?

If you can't answer this question, it's worth considering this product. Easily fitted to an O<sub>2</sub> therapy door the Handi+ displays the concentration level on an LCD display. It can also be used in line with your anaesthetic machine, ideal for verifying output from O<sub>2</sub> concentrators.

It is simple to use, with one-touch calibration, a new over-mold, calibration reminder, and auto-off, combined with a 2-year warranty to make it the best of its class.

#### Features

- New and improved one-touch calibration button
- Durable over-mold for extra protection
- Small, convenient, and portable
- Auto-off
- Lanyard Attachment
- Protective Wrap

### ADDITIONAL INFORMATION

<b>Animal</b>	Mixed
<b>SKU</b>	129-555
<b>Colour/Style</b>	Blue
<b>Monitor Parameters</b>	Oxygen

<b>Screen Type</b>	Digital Display
<b>Power Source</b>	Battery
<b>Optional Modules</b>	No
<b>SD Card</b>	No

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.  
**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!  
**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.  
**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights