

**VET BP DOPPLER KIT**

DOP-143

**PRODUCT DESCRIPTION****Burtons Vet BP Doppler cased kit complete with 8 MHz probe sphygmomanometer & 4 cuffs + gel**

The Doppler Vet BP is specifically designed for a reliable and accurate blood pressure measurement on small animals. easy to set up and use the Vet BP is allows for easy flow detection in the artery.

The probe can be fixed on the patient paw or tail during anaesthesia.

**Features**

- Sensitive, wide beam, focused 8 MHz probe for better flow detection
- Built in rechargeable battery
- Hand-held design
- Optional pencil probe for exotics
- Optional transesophageal probe for cardiac monitoring

**Set includes:**

1x Doppler

1x transversal linear probe 8 MHz

1x battery charger

1x sphygmomanometer

4x cuffs in different sizes

1x headset

1x tube of gel

### Specifications:

- Ultrasonic frequency : 8 MHz
- Audio output : 200mW
- Frequency response : 300Hz - 6kHz
- Battery : 3 x 1.2V, 800mAh €ç Operating time : +/- 6 hours
- Charging time : 15 hours
- Dimensions : 175 x 80 x 25mm
- Weight incl. battery : 300 g
- Headset acoustic output : mini jack

### ADDITIONAL INFORMATION

<b>Animal</b>	Small Animal
<b>SKU</b>	DOP-143
<b>Overall Product Dimensions (mm)</b>	175 x 80 x 25mm
<b>Colour/Style</b>	Black
<b>Brand</b>	Burtons
<b>Monitor Parameters</b>	NIBP
<b>Screen Type</b>	N/A
<b>Power Source</b>	Rechargeable
<b>Optional Modules</b>	No
<b>SD Card</b>	No

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.  
**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!  
**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.  
**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

