

## OESOPHAGEAL PROBE FOR VET BP DOPPLER

DOP-144



### PRODUCT DESCRIPTION

## Add cardiac monitoring to your Burtons Vet BP Doppler

Additional Oesophageal probe suitable for Burtons Vet BP Doppler

### ADDITIONAL INFORMATION

|                     |                |
|---------------------|----------------|
| <b>Animal</b>       | Small Animal   |
| <b>SKU</b>          | DOP-144        |
| <b>Colour/Style</b> | Black          |
| <b>Brand</b>        | Burtons        |
| <b>Type</b>         | Probe          |
| <b>Monitor Fit</b>  | Vet BP Doppler |

## Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights