


## CO<sub>2</sub> MONITOR EXHAUST KIT

V1175



### PRODUCT DESCRIPTION

## Enable sample gases to be evacuated into your scavenge system

- CO<sub>2</sub> Monitor Exhaust Kit
- Complete with Scavenging tube (Adaptor V1175-003) and Scavenge Connector 22mm M/F with port (V1175-003)
- For monitors: V8200, V8401, V8402, V90041, V90043 V9400 Advisor Monitor & All IMEC's

### ADDITIONAL INFORMATION

<b>Animal</b>	N/A
<b>SKU</b>	V1175
<b>Brand</b>	Burtons
<b>Type</b>	Other
<b>Monitor Fit</b>	Cardell 9401, Surgi Vet Advisor

## Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights