

**GAS MANIFOLD FILTER 10 PK**

V8208



Category: Patient Monitoring &gt; Monitoring Accessories

Price: £101.87

Date Printed: 28/03/24

**PRODUCT DESCRIPTION****SurgiVet CO<sup>2</sup> Gas Manifold Filters**

- Suitable for use with V8400 Series Hand Held Capnographs.
- Also useful as a filter in between the sample line and moisture trap on the Advisor monitor
- Pack of 10 filters

**ADDITIONAL INFORMATION**

|                     |                   |
|---------------------|-------------------|
| <b>Animal</b>       | Small Animal      |
| <b>SKU</b>          | V8208             |
| <b>Colour/Style</b> | Clear             |
| <b>Brand</b>        | SurgiVet          |
| <b>Type</b>         | Other             |
| <b>Monitor Fit</b>  | Surgi Vet Advisor |

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights